
api®

BIOMÉRIEUX

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

API 20 E

+ - + - + -

Tests $\leq$ 3
(y compris / including
test GLU)

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






操作步骤
Procedure

References

Based on numerical identification concept, API® has made microbial identification easy, rapid and reliable.

Unrivalled for over thirty years in terms of performance and the range available, it is now considered as the reference :

- Used to evaluate any new identification system arriving on the market.
- Used to teach microbiological identification in most higher education centres.

Reliability

Above all API offers reliability, combining a simple and standardized method with extensive and extremely robust databases for all microorganisms. For this reason it is easily inserted into a quality procedure.

Savings

The lack of specific investment, very short training period and simplified order management due to its very long shelf life all make API range a very economical routine identification method.

api®


BIOMÉRIEUX

Gram⁻ bacilli

Orientation tests

Other Gram⁻ bacilli

Orientation tests

Gram⁻ bacilli

Gram⁻ cocci in pairs (*Neisseria*)
Small Gram⁻ bacilli (*Haemophilus*)

Growth with CO₂ atmosphere

api NH

Neisseria haemophilus

Curved bacilli,
or in spiral

Growth with microaerophilic atmosphere

Oxydase ⁺

api CAMPY

Campylobacter


BIOMÉRIEUX

Gram⁺ bacilli

Orientation tests


BIOMÉRIEUX

Gram⁺ cocci

Orientation tests

Yeasts

Orientation tests

Microscopic observation :
YEASTS

Assimilation tests

Lapi 20C AUX

Candida
Cryptococcus
Geotrichum
Kloeckera
Pichia
Rhodotorula
Trichosporum
Saccharomyces

Sugar acidification or
enzymatic tests

Lapi CANDIDA

Candida
Cryptococcus
Saccharomyces
Trichosporum

Reagents

	API 20 E* (a)		RapiD 20E		API 10 S		API 20 NE		API Staph		API 20 Strep		API Coryne		API Listeria		API Candida		API 20 C AUX		API 20 A*		API Campy		API NH		API ZYM		API 50 CH		API 50 CHB/E		API 20 E (a)		API 50 CH		
	20100	20700	10100	20050	20500	20600	20900	10300	10500	20210	20300	20800	10400	25200	50300	50430	20100	50300	50410																		
ADDITIONAL REAGENTS																																					
Zn (2 x 10 g)	x			x																																	
TDA (2 ampoules)	1/8		1/4																																		
VP 1 VP 2 (2 x 2 ampoules)	1/8	1/8			1/8	1/8																															
NIT 1 NIT 2 (2 x 2 ampoules)	1/8		1/4	1/8	1/8		1/20							1/20																							
ZYM A (2 ampoules) ZYM B (2 ampoules)					1/8 1/8	1/8 1/8	1/20 1/20																			1,5 1,5											
NIN (1 ampoule)						1/2															1/2																
BCP (1 ampoule)																				1/2																	
EHR (1 ampoule)																				1/2																	
XYL (2 ampoules)																				1/8																	
JAMES (2 ampoules)	1/8	1/8	1/4	1/8																																	
FB (2 ampoules)																					1/4																
PYZ (1 ampoule)							1/10																														
SUSPENSION MEDIA																																					
NaCl 0.85% medium (5ml)																																					
NaCl 0.85% medium (3ml)																																					
NaCl 0.85% medium (2ml)			1/4		1/4																																
Suspension medium (5ml)	1/4		1/2																																		
Suspension medium (2ml)						1/4										1/4																					
ADDITIONAL PRODUCTS																																					
Oxidase Reagent	x	x	x	x																																	
API OF medium	x	x	x	x																																	
API M medium	x	x	x	x						x																											
Paraffin oil	x	x	x	x	x	x	x	x			x					x			x	x	x	x	x	x													
Rack of 12 ampoules	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
Sterile PSlipettes (5 ml)	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
Sterile swabs							x	x								x			x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		
Mc Farland Standard		x		x	x	x			x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	x		

*figures = recommended quantity of additional reagent(s) for each strip kit ordered

(eg. TDA = 1/8 <=> 1/8 of kit ref. 70402 (2 TDA ampoules) is required to use one kit of ref. 20100 (25 API 20 E strips))

x = additional product required

(a) : for the API 20 E strip, use the additional reagent kit ref. 20120 (7 ampoules)

BIOMÉRIEUX

**Simplified
methodologies**

	Ref.	Microorganisms	Suspension	Mc Farland	Transfert	Medium	Incubation	Atmosphere
API 20 E *	20100 20160	<i>Enterobacteriaceae</i> and other non-fastidious Gram negative bacilli	5 ml NaCl 0.85% medium or 5 ml Suspension medium	1 colony	NA	NA	37°C 18-24hrs/48hrs	Aerobic
RapiD 20E	20701	<i>Enterobacteriaceae</i>	2 ml NaCl 0.85% medium	0.5 McF	NA	NA	37°C 4hrs	Aerobic
API 10 S	10100	<i>Enterobacteriaceae</i> and other non-fastidious Gram negative bacilli	5 ml NaCl 0.85% medium or 5 ml Suspension medium	1 colony	NA	NA	37°C 18-24hrs	Aerobic
API 20 NE	20050	Non- <i>Enterobacteriaceae</i> and non-fastidious Gram negative bacilli	2 ml NaCl 0.85% medium	0.5 McF	200 µl	AUX medium	30°C 24-48hrs	Aerobic
API Staph	20500	<i>Staphylococcus</i> , <i>Kocuria</i> <i>Micrococcus</i> genera	API Staph medium	0.5 McF	NA	NA	37°C 18-24hrs	Aerobic
API 20 Strep	20600	Streptococci	2 ml Suspension medium	4 McF	500 µl	API GP medium	37°C 4-24hrs	Aerobic
API Coryne	20900	Coryneform bacteria	3 ml Suspension medium	6 McF	500 µl	API GP medium	37°C 24hrs	Aerobic
API Listeria	10300	<i>Listeria</i>	2 ml Suspension medium	1 McF	NA	NA	37°C 18-24hrs	Aerobic
API Candida	10500	Yeasts	2 ml NaCl 0.85%	3 McF	NA	NA	37°C 18-24hrs	Aerobic
API 20 C AUX	20210	Yeasts	2 ml NaCl 0.85% 2 ml Suspension medium	2 McF	100 µl	API C medium	30°C 48-72hrs	Aerobic
API 20 A *	20300	Anaerobes	API 20 A medium	3 McF	NA	NA	37°C 24hrs	Anaerobic
API Campy	20800	<i>Campylobacter</i>	3 ml NaCl 0.85% medium	6 McF	150 µl	AUX medium	37°C 24-48hrs	Aerobic/CO ₂
API NH	10400	<i>Neisseria</i> , <i>Haemophilus</i> and <i>Branhamella catarrhalis</i>	2 ml NaCl 0.85%	4 McF	NA	NA	37°C 2hrs	Aerobic
API ZYM	25200	Semi-quantitative enzyme activity tests	2 ml NaCl 0.85% or 2 ml Suspension medium	5-6 McF	NA	NA	37°C 4hrs	Aerobic
API 50 CH API 50 CHB/E	50300 50430	<i>Bacillus</i>	5 ml NaCl 0.85% + 50 CHB/E medium	2 McF 2 McF	NA	NA	30°C 24-48hrs 55°C 3/6-24hrs	Aerobic
API 50 CH API 50 CHB/E	50300 50430	<i>Enterobacteriaceae</i>	NaCl 0.85% 5 ml + 50 CHB/E medium	0.5 McF 0.5 McF	NA	NA	37°C 24-48hrs	Aerobic
API 50 CH API 50 CHL	50300 50410	<i>Lactobacillus</i>	50 CHL medium	2 McF	NA	NA	30/37°C 48hrs	Aerobic

BIOMÉRIEUX

anytime, anywhere

Simplified interpretation for every lab

By miniaturizing and standardizing conventional techniques the API[®] range has represented a true revolution in the field of identification.

Committed to remaining at the forefront of innovation, bioMérieux once again gives a new dimension to this range by now making identification available to everyone, at anytime, and anywhere.

APIWEB[™], a tool which encompasses the technologies of today and the future, is the essential complement to manual identification.

User-friendly

- Internet platform
- Easy-to-use tool, integrated intuitive interface using Internet navigation tools
- Available 7 days a week, 24 hours a day
- Reading assistance for each identification strip

High-performance

APIWEB contains all the regularly updated API/ID32 databases. Containing more information than API Analytical Profile Indexes, APIWEB interprets a great number of profiles enabling the identification of 600 species of bacteria and yeasts

- it is a reliable, state-of-the-art interpretation tool.

Each identification is associated with a complete report :

- species proposed
- comments reflecting the reliability of the identification, selected after the calculation of two indices :
 - the identification percentage (probability of species identification)
 - and the typicity index (typical character of the profile studied)
- complete biochemical profile
- additional information and tests

Security

- Secure site, personalized access codes
- Data encoding ensuring the integrity and confidentiality of your results
- Highly reactive tool maintenance

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