


helica Ochratoxin A ELISA

For Wine And Grape Products



[Catalog No. 981OCH01W-48, 981OCH01W-96, 981OCH01ALC]

Product Highlights

- Visual reading of 0, 1, and 2 ppb
- Limit of detection < 1ppb
- Easy extraction
- Cost-effective and user-friendly tests
- Results in less than 30 minutes
- Break-away wells allow low-volume testing
- Standards and conjugate preserved with proprietary stabilizers
- One year stability and reproducibility of results
- Easy to follow instructions
- No confirmation by HPLC required

Mycotoxin ELISA KITS

Catalog No.	Product
961AFLM01M-96	Aflatoxin M1 for milk and dairy products
941BAFL01B-96	Aflatoxin B1 for grains, cereal, nuts, & animal feed
941AFL01M-96	Total Aflatoxins for grains, cereal, nuts, & animal feed
981AFL01LM-96	Total Aflatoxins Low Matrix for grains, cereal, silage, nuts, spices, & animal feed
991AFLM01U-96	Aflatoxin M1 For Urine
941DON01M-96	Deoxynivalenol (DON) for grains, cereal, & animal feed
951FUM01C-96	Fumonisin for maize (corn)
951ZEA01G-96	NEW Zearalenone for corn, barley, oats, wheat, rice, sorghum
941OCH01M-96	Ochratoxin A For grains, cereals, animal feed- EU Standards
981OCH01W-96	Ochratoxin A Qualitative for wine, grape juice & grape must
981OCH01ALC-96	Ochratoxin A Quantitative for alcoholic beverages, wine, grape juice, grape must, & beer
991OCH01MS-96	Ochratoxin A for human and animal milk, serum, & plasma
961OCH01COF-96	Ochratoxin A for coffee, cocoa butter and spices



Helica Biosystems, Inc.
1527 W. Alton Ave
Santa Ana, CA 92704, USA

P: 714-578-7830

F: 714-578-7831

www.helica.com

info@helica.com